

# Bailey

SURGICAL INSTRUMENT  
SPECIALISTS

.....

## EXPORT PODIATRY CATALOGUE EDITION 3



[BaileyInstruments.co.uk](http://BaileyInstruments.co.uk)



MD642415

**Supporting Podiatry for over 30 years with high quality products and devices**

**We are Global Specialists in Podiatry Instruments and Devices.**

Whether it is reusable instruments, single use or diagnostic devices for diabetes you can rely on us for innovation, quality, precision, performance and value for money. With our speedy and reliable logistics from our Global Supply Centre based in Manchester UK we can supply any market.



Contact Gemma Bailey, our Export Director on **+44 7545 915 578**  
Email: [export@baileyinstruments.co.uk](mailto:export@baileyinstruments.co.uk)

To see our full range of products, please visit our website: [www.baileyinstruments.co.uk](http://www.baileyinstruments.co.uk)

Serving the Podiatry profession for over 30 years, Bailey Instruments Ltd is a leading contract supplier to the NHS throughout the UK.



NHS Contract Supplier for over 30 years

Global Supplier to over 30 countries



Account Number:

---

---

---

---

---

---

---

---

---

---

**Hardenburg Range** 4 - 6

**Bailey Reusable**

- Nippers 7 - 10
- Blacks Files 10
- Scalpels and Blade Safety 11
- Swann Morton Blades 12
- Probes and Files 13
- Spatulas and Burs 14 - 15
- Scissors/ Nail Surgery / Nippers 16 - 18
- Skin Surgery 19

**Susol Single Use**

- Podiatry Packs 20
- Podiatry Instruments 21 - 23

**Diagnostics**

- Bailey Monofilaments 24
- Neuropathy 25
- Patella Hammers & Tuning Forks 26
- Infrared Thermometer 27

**Orthotics**

- Bailey Insoles 28 - 29

**Miscellaneous Instruments**

- Scissors, Forceps & Reach Your Feet 30
- How To Trade Internationally 31

REUSABLE - HARDENBURG RANGE






**HARDENBURG**  
S O L I N G E N

Manufactured from high grade German steel, our Hardenburg range is well recognised for its quality and performance. New additions to the range include diabetic and cantilever nippers.

REUSABLE - HARDENBURG RANGE

**Hardenburg General Nipper**  
CURVED, SCREW JOINT




Code	Size
CH899B	14cm

**Hardenburg Thwaites Nipper**



Code	Size
CH890	13cm

**Hardenburg General Nipper**  
CURVED, SQUARE HANDLES, SCREW JOINT



Code	Size
CH892	14cm

**Hardenburg General Nipper**  
CURVED, BOX JOINT, CATCH



Code	Size
CH893	14cm

**Hardenburg Honeback Nipper**  
STRAIGHT, BOX JOINT



Code	Size
CH895	14cm

**Hardenburg Badger Handle**



Code
CH891

**Hardenburg General Nipper**  
STRAIGHT, BOX JOINT, CATCH



Code	Size
CH894	14cm

**Hardenburg General Nipper**  
STRAIGHT, SCREW JOINT



Code	Size
CH899A	14cm

**Hardenburg Ingrown Nipper**  
BOX JOINT



Code	Size
CH896	11cm

**Hardenburg Ingrown Nipper**  
BOX JOINT




Code	Size
CH897	13cm

REUSABLE - HARDENBURG RANGE

**Hardenburg Nipper Ingrown**

BOX JOINT




Code	Size
CH898	15cm

**Hardenburg Diabetic Pedicure Nipper**



Code	Size
CH902	14cm

**Hardenburg Diabetic Corner Nipper**



Code	Size
CH901	13cm

**Hardenburg Cantilever Convex Jaw**



Code	Size
CH900	16cm



**Bailey**

SURGICAL INSTRUMENT SPECIALISTS

**REUSABLE INSTRUMENTS**

We have been manufacturing podiatry instruments for over 30 years. Our range of reusable instruments are used the world over and recognised for their quality and durability.



**General Nipper**

CURVED, SCREW JOINT, BARREL SPRING



Code	Size
CH01	14cm

**Ingrown Nipper**

STRAIGHT, BOX JOINT



Code	Size
CH02	12.5cm
CH02A	15cm

**General Nipper**


CURVED, BOX JOINT, CATCH



Code	Size
CH03	14cm

**General Nipper**

STRAIGHT, BOX JOINT



Code	Size
CH30	14cm



BAILEYS REUSABLE - NIPPERS

**General Nipper**


MOON END CUT



Code
CH112

**Bevelled General Nipper**

STRAIGHT, BOX JOINT, CATCH



Code	Size
CH32A	14cm

**General Nipper**

CURVED, SCREW JOINT, CATCH



Code	Size
CH35A	14cm

**Ultra-Lite Nipper**


STRAIGHT / CURVED



Code	Size	Type
CH90	12.5cm	STR
CH90A	15cm	STR
CH91	12.5cm	CVD
CH91A	15cm	CVD

**Bevelled General Nipper**


CURVED, BOX JOINT



Code	Size
CH32B	14cm

**General Nipper**

STRAIGHT, SCREW JOINT



Code	Size
CH34	14cm

**Cantilever Nipper**

CURVED, CATCH



Code	Size	Catch
CH23A	15cm	Yes
CH23B	15cm	No

**Knurled Handle Cantilever Nipper**

STRAIGHT / CURVED



Code	Size	Type
CH23S	15cm	STR
CH23	12.5cm	CVD

**General Nipper**


STRAIGHT, SCREW JOINT, CATCH



Code	Size
CH34A	14cm

**General Nipper**

CURVED, SCREW JOINT



Code	Size
CH35	14cm

**Cantilever Nipper**

CURVED, DOUBLE SPRING



Code	Size
CH23C	15cm

**Cantilever Nipper**

CURVED, FINE HEAD, CATCH



Code	Size
CH23D	15cm

**BAILEYS REUSABLE - NIPPERS & BLACKS FILES**


**BAILEYS REUSABLE - SCALPELS AND BLADE SAFETY**

**Liston Cantilever Nipper**  
ANGLED, STRAIGHT



Code	Size
B117	14cm

**Cantilever Nipper**  
CURVED, SILICONE GRIP



Code	Size
CH23G	15cm

**Blacks File**  
S/E FINE CUT



Code	Size
CH07	13cm

**Blacks File**  
S/E STANDARD CUT



Code	Size
CH07A	13cm

**Blacks File**  
D/E FINE / STANDARD



Code	Size
CH07B	14cm

**Blacks File**  
S/E MICRO HEAD



Code	Size
CH07C	14cm

**Blacks File**  
D/E WITH PROBE



Code	Size
CH07E	14cm

**Blacks File**  
D/E ANGLED



Code	Size
CH07F	15cm

**Blacks File**  
D/E CURVED



Code	Size
CH07G	14cm

**Blacks File**  
S/E LARGE HEAD



Code	Size
CH07H	14cm

**Badger Handle**  
KNURLED



Code
CH105

**Badger Handle**  
HEXAGONAL



Code
CH105A

**Swann-Morton® Scalpel Handle**  
NO. 9




Code
CH209

**Swann-Morton® Scalpel Handle**  
NO. 3 GRAD



Code
CH56G

**Scalpel Handle**  
BAILEY NO.3



Code
CH56BG

**Swann-Morton® Scalpel Handle**  
NO. 4 GRAD



Code
CH64G

**Forceps**  
SCALPEL BLADE REMOVING



Code	Size
CH18	16cm

**Clicksmart BladeFLASK**

Clicksmart BladeFLASK provides a simply safer way to remove and contain used scalpel blades.



Code
SM-QYUKSM

**Swann-Morton® Disposable Handle with Blade**

Code	Description
SS04-6601	No.3 Handle with No.10 Blade
SS04-6603	No.3 Handle with No.11 Blade
SS04-6605	No.3 Handle with No.15 Blade



**Swann-Morton® Blade Box**  
HOLDS 100 BLADES

Code
CN8B



Swann-Morton® **BLADES**

Swann-Morton® **NO.10**



Code	Unit
SM10	100 - Sterile

Swann-Morton® **NO.10A**



Code	Unit
SM10A	100 - Sterile

Swann-Morton® **NO.11**



Code	Unit
SM11	100 - Sterile

Swann-Morton® **NO.E/11**



Code	Unit
SME/11	100 - Sterile

Swann-Morton® **NO.E/11S**



Code	Unit
SME/11S	100 - Sterile

Swann-Morton® **NO.15**



Code	Unit
SM15	100 - Sterile

Swann-Morton® **NO.15T**



Code	Unit
SM15T	100 - Sterile

Swann-Morton® **FINE NO.61**



Code	Unit
SM61	25 - Sterile

Swann-Morton® **FINE NO.62**



Code	Unit
SM62	25 - Sterile

Swann-Morton® **FINE NO.64**



Code	Unit
SM64	25 - Sterile

Swann-Morton® **FINE NO.65**



Code	Unit
SM65	25 - Sterile


Swann-Morton® **FINE NO.67**



Code	Unit
SM67	25 - Sterile


**BAILEYS REUSABLE - PROBES & FILES**

**Probe Swan Neck**  
SPOON/HEAD D/E




Code
CH102M

**Probe Round Head**  
WITH HOLE




Code
CH206

**Probe Round Head**  
SWAN NECK



Code
CH206M

**Probe Swan Neck**  
FLAT



Code
CH25

**Probe Swan Neck**  
FLAT HEAD



Code
CH26
CH26L

**Foot Dresser / File**



Code	Size
CH140	20cm

**Foot Dresser / File**



Code	Size
CH140A	15cm

**Foot Dresser / File**  
NARROW




Code	Size
CH140B	15cm

**Foot Dresser / File**  
PLASTIC PADDLE



Code
CH143

**Foot Rasp**



Code
CH39

**Tube Gauze Applicator**



Code	Size
CH100	00

**Tube Gauze Applicator**



Code	Size
CH101	01

## BAILEYS REUSABLE - SPATULAS AND BURS

### Packer / Spatula KNURLED HANDLE



Code  
CH207

### Spatula with Knurled Grip



Code	Length
CH27B	12.5cm

### Spatula 12.5CM



Code  
CH27

### Spatula with Hexagonal Grip



Code	Length
CH27A	16.5cm

## Reusable Burs

### Bur BUD 1



Code  
DB2

### Bur BUD 2




Code  
DB3

### Bur PEAR 1



Code  
DB4

### Bur DIAMOND 15BHP



Code  
DB5

### Bur BUD 3



Code  
DB17

### Bur PEAR 2



Code  
DB18

### Bur FLAME



Code  
DB19

### Bur CYLINDER 1



Code  
DB1

### Bur 15C



Code  
DB7

### Bur 15HP



Code  
DB6

### Bur 102HP



Code  
DB11

## BAILEYS REUSABLE - BURS

### Bur 101HP



Code  
DB10

### Bur 103HP



Code  
DB12

### Bur 104HP



Code  
DB13

### Bur 105HP



Code  
DB14

### Bur 106HP



Code  
DB15

### Bur 1HP



Code  
DB8

### Bur 2HP



Code  
DB9

### Bur TUNGSTEN F28




Code  
DBF28

### Bur TUNGSTEN F29



Code  
DBF29

### Bur TUNGSTEN F30




Code  
DBF30

### Bur CERAMIC




Code  
DBC31

### Mandrel Clip




Code  
DB16

### Autoclave Bur Box OVAL




Code  
CH99

### Autoclave Bur Box SLIDING



Code  
CH99A

### Autoclave Bur Box CASKET



Code  
CH99B

### Sanding Discs

Code	Size
CN-SD19C	19mm Coarse
CN-SD19M	19mm Medium
CN-SD22C	22mm Coarse
CN-SD22M	22mm Medium





**BAILEYS REUSABLE - SCISSORS / NAIL SURGERY**

**Dressing Scissors**  
BL/SH STR



Code	Size
CH10	14cm
CH10L*	14cm


\*L = Left Handed

**Scissors Dressing**  
BL/SH STR 15CM



Code	Size
CH69	15cm
CH69L	18cm

**Dressing Scissors**  
BL/SH STR TUNGSTEN CARBIDE



Code	Size
CH69TC	15cm

**Susol Iris Scissors**  
FINE POINT STR 11.5CM



Code	Type
CH13	Straight
CH14	Curved

**Mayo Scissors**  
STRAIGHT



Code	Size
SD7-0300	14cm
SD7-0301	16.5cm
SD7-0302	19cm


**Mayo Scissors**  
CURVED



Code	Size
SD7-0305	14cm
SD7-0306	16.5cm
SD7-0307	19cm

**BAILEYS REUSABLE - NIPPERS / NAIL SURGERY**

**Thwaites Nipper**  
ENGLISH PATTERN



Code	Size
CH21A	14cm

**Thwaites Nipper**  
PISTOL GRIP



Code	Size
CH21B	12.5cm

**Thwaites Nipper**  
LISTON




Code	Size
CH21C	15cm

**Thwaites Nipper**  
ROUNDED




Code	Size
CH21D	13.5cm

**Nail Elevator**  
LOCKE



Code
CH208

**Nail Elevator**  
FRANCIS



Code
CH203

**Nail Splitter**  
FRANCIS



Code
CH204

**Trocar Handle**  
LOCKE ELEVATOR




Code
CH208A

**BAILEYS REUSABLE - NAIL SURGERY**

**Artery Forceps**

HALSTEAD MOSQ STRAIGHT

Code	Size
SD9-0625	12.5cm



**Artery Forceps**

HALSTEAD MOSQ CURVED

Code	Size
SD9-0627	12.5cm

**Artery Forceps**

SPENCER WELLS STRAIGHT

Code	Size
SD9-0835	12.5cm



**Artery Forceps**

SPENCER WELLS CURVED

Code	Size
SD9-0845	12.5cm

**Ramplays Sponge Holding Forceps**



Code	Size
SD16-1481	15cm
SD16-1482	18cm
SD16-1483	24cm

**Needle Holder**


MAYO HEGAR



Code	Size
SD18-1770	12.5cm

**Adson Dissecting Forceps**

SERRATED



Code	Type
SD11-1020	Serrated
SD11-1021	1:2 Teeth*

\* Not pictured


**Scissors**

LISTER BANDAGE



Code	Size
SD7-0291	14cm
SD7-0292	18.5cm

**McDonalds Dissector**



Code	Size
SD49-3505	19cm

**BAILEYS REUSABLE - SKIN SURGERY**

**McIndoe Dissecting Forceps**



Code	Size
SD11-1264	15cm

**Kilner Catspaw Retractor**

D/E



Code	Size
SD25-2262	15cm

**Dissecting Forceps**

TOE



Code	Size
SD11-0951	15cm

**Dissecting Forceps Gillies**

1:2 TEETH



Code	Size
SD11-1180	15cm

**Allis Tissue Forceps**

4:5 TEETH



Code	Size
SD50-3551	15cm

## SUSOL SINGLE USE - PODIATRY PACKS



### SINGLE USE PODIATRY TOOLS

This range of instruments was designed to reduce patient risk, be cost effective and increase productivity. Our unique and patented nipper gives you a sharp cut every time, while improving comfort through the shock absorbent silicone handles.



### Susol Podiatry Set

- Nipper 13cm CVD
- Blacks File 13cm STD
- Scalpel Handle No.3
- Dressing Scissors 11.5cm
- Emery Board 18cm



Code	Qty
<b>BSDP-01A</b>	<b>100</b>
<b>BSDP-01P</b>	<b>10</b>

### Susol Nail Care Set

- Nipper 13cm CVD
- Blacks File 13cm STD
- Emery Board 18cm



Code	Qty
<b>BSDP-02</b>	<b>100</b>
<b>BSDP-02P</b>	<b>10</b>

### Susol Basic Care Set

- Nipper 13cm CVD
- Emery Board 18cm



Code	Qty
<b>BSDP-04</b>	<b>100</b>
<b>BSDP-04P</b>	<b>10</b>

### Susol Nail Surgery (PNA) Set

- Thwaites Nipper
- Beaver Handle
- Locke Elevator
- Mosquito Forceps
- Blacks File
- Scalpel Handle No.3
- Dressing Scissor
- Ribbon Tourniquet



Code	Qty
<b>BSDP-03</b>	<b>10</b>

### Wound Debridement Kit

- Eliminates the need to search for missing procedural components or the waste of extra unused items
- Helps to contain costs and promote standardised procedures
- Encourages consistent product usage and compliance with clinical protocol



Code	Qty
<b>SSP-130</b>	<b>10</b>

## SUSOL SINGLE USE - PODIATRY INSTRUMENTS

### Susol General Nipper

CVD 13CM

- Supercut blades
- Concealed spring
- Shock absorbing handles



Code	Qty
<b>SS08-0500</b>	<b>10</b>

### Susol Ingrown Nail Nipper

14CM

- Supercut blades
- Concealed spring
- Shock absorbing handles



Code	Qty
<b>SS08-0537</b>	<b>10</b>

### Susol Thwaites Nipper


14CM



Code	Qty
<b>SS08-0539</b>	<b>10</b>

### Susol Artery Forceps

HALSTEAD MOSQ 12.5CM



Code	Qty	Type
<b>SS09-0625</b>	<b>10</b>	<b>STR</b>
<b>SS09-0627</b>	<b>10</b>	<b>CVD</b>

### Susol Dressing Scissor

PLASTIC HANDLE 11.5CM



Code	Qty
<b>BSD-01</b>	<b>10</b>

### Susol Dressing Scissor

B/B STR 13CM




Code	Qty
<b>SS07-0234</b>	<b>10</b>

## SUSOL SINGLE USE - PODIATRY INSTRUMENTS

### Susol Iris Scissors

FINE POINT STR 11.5CM



Code	Qty
SS07-0282	10

### Susol Lister Bandage Scissors



Code	Size	Qty
SS07-0289	11.5cm	10
SS07-0290	14cm	10
SS07-0291	18.5cm	10

### Susol Digital Tourniquets



Code	Size	Qty
S-TQS	Small	50
S-TQM	Medium	50
S-TQL	Large	50

### Susol Tube Gauze Applicator



Code	Size	Qty
SS08-0561	Small	50
SS08-0562	Large	50

## Susol Disposable Burs

### Susol Bur

BUD 1



Code	Qty
S-DB2	100
S-DB2P	50

### Susol Bur

BUD 2



Code	Qty
S-DB3	100
S-DB3P	50

### Susol Bur

BUD 3



Code	Qty
S-DB17	100
S-DB17P	50

### Susol Bur

PEAR 1



Code	Qty
S-DB4P	100
S-DB3P	50

### Susol Bur

PEAR 2



Code	Qty
S-DB18	100
S-DB18P	50

### Susol Bur

IHP



Code	Qty
S-DB8	100
S-DB8B	50

### Susol Bur

15C



Code	Qty
S-DB7	100
S-DB7P	50

### Susol Wound Depth Gauge



Code	Qty
S-DG1	100

### Susol Emery Board

18CM



Code	Qty
CH145	100

### Susol Blacks File SE



Code	Cut	Qty
SS08-0552	Standard	10
SS08-0553	Coarse	10

### Susol Blacks File

MICRO HEAD WITH ROUND PROBE



Code	Qty
SS08-0557	10

### Susol Scalpel Handle



Code	Size	Qty
SS04-0100	NO.3	10
SS04-0106	NO.9	10

### Susol Adsons Dissecting Forceps



Code	Size	Qty	Type
SS11-1020	13cm	10	Toothed
SS11-1021	13cm	10	Non-Toothed

### Susol McKays Elevator



Code	Size	Qty
SS08-0567	16.5cm	10

### Susol Locke Elevator



Code	Size	Qty
SS08-0548	14cm	10

### susol House Curette

DOUBLE ENDED, WITH HOLE



Code	Size	Qty
SS62-4300	180mm	10

# Bailey

## THE BAILEY RETRACTABLE ACCURATE MONOFILAMENT

Now regarded as a cornerstone in the screening and assessment of the diabetic foot. The Bailey Retractable Monofilament is used widely throughout the world.



## The Cornerstone of Screening and Assessment

The most important feature of our monofilament is its accuracy. A patient who is unable to feel a 10g monofilament on their foot is said to have lost protective sensation. This increases the patients' risk of developing a foot ulcer so an accurate monofilament is essential.

The use of an accurate monofilament is recommended in the NG19 Guideline for Diabetic Foot Problems (National Institute for Clinical Excellence), SIGN 116 management of diabetes (Scottish Intercollegiate Guidelines Network) and the ADA (American Diabetes Association).

The Bailey Retractable Monofilament is accurate for 10,000 'bounces' equating to 1000 patients. It is the most economical tool for use in screening and assessment. For more information, please visit: [www.AccurateMonofilament.com](http://www.AccurateMonofilament.com)

### Monofilament

10G

Code  
CH537



### Monofilament

75G

Code  
CH534



## Neurothesiometer

Code	Description
NEU1501	Unit with Battery
NEU1520	Spare Battery
NEU1500	Calibration Service

- Battery operated diagnostic instrument for the assessment of vibration thresholds on the surface of the body.
- Housed in a strong, compact carrying case.
- The Neurothesiometer has an LED digital display.

Code	Description
NEU1522	UK Charger (3 pin)
NEU1523	USA Charger
NEU1524	Continental Charger (2 pin)
NEU1525	South Africa Charger
NEU1526	Australia Charger



## Two Point Discriminator

Used for 2 point discrimination test to assess tactile perception.

Code  
CH71



## Tip Therm

This unique tool measures temperature sensitivity of the skin. Made of polymer and metal alloys. The polymer end feels warmer, the alloy end feels cooler due to the thermal conductivity properties of these materials.

Code  
CH536



## DIAGNOSTICS - PATELLA HAMMERS, TUNING FORKS

### Patella Hammer

BRUSH & PIN



Code  
CH80

### Patella Hammer

NYLON



Code  
CH81

### Patella Hammer

TELESCOPIC



Code  
CH82

### Bailey Patella Hammer

WITH PIN



Code  
CH83

### Tuning Fork

Used for checking vibration sensation in examination of peripheral nervous system



Code	Frequency
CH128	128Hz
CH256	256Hz

### Tuning Fork

RYDEL SEIFFER

Vibrates at 64Hz and ideal for detecting impaired vibration sensitivity



Code  
CH257

## DIAGNOSTICS - INFRARED THERMOMETER

### Infrared Thermometer

Ideal for comparative foot temperature measurement. A great tool in helping to establish the risk of foot ulceration in the diabetic patient.

- Easy to use
- Clear display
- Reading memory feature
- Up to 1,000 measurements on one battery set
- Easily cleaned
- Accurate - Max variance of +/-0.3 centigrade
- Batteries and manual included

#### SUGGESTED USES:

- For temperature screening for pre-entry to buildings
- To indicate infection sites arounds wounds
- Charcot foot screening/examination

Code  
D/A201



NEW FROM Bailey




We've introduced two new ranges to our orthotic offerings. Our Versatile insoles offer a softer, kinder material for improved patient comfort and compliance.

### Bailey Versatile

The material offers a softer, kinder control leading to improved patient compliance. Infinitely adjustable, the insoles can be trimmed to fit and can accept all additions, such as wedges and heel raises. Ideal for at risk feet and perfect for support in high impact sports.

Code	Sizing		
VER-A	Small	UK 3-5	EU 36-38
VER-B	Medium	UK 6-8	EU 40-42
VER-C	Large	UK 9-11	EU 43-46



### Bailey Versatile 3/4


The three quarter length device is designed for quick and easy fitting into most footwear. A low shore rating provides more comfort for the user and increases the likelihood of compliance and positive outcome.

Code	Sizing		
VERC-A	Size 3-5	EU 36-38	EU 36-38
VERC-B	Size 6-8	EU 40-42	EU 40-42
VERC-C	Size 9-11	EU 43-46	EU 43-46



### Pedum Flex

The Pedum Flex mirrors the Firm but with a lower shore rating of 70. Heel raises and medial heel wedges can be used facilitating rapid clinical adjustment.



Code	Size
FLEX	From shoe size Jnr 12 - adult 11'

### Pedum 5 Degree Lateral

Designed for Osteoarthritis of the knee



Size
From shoe sizes 3-11

### Pedum Firm

With a shore rating of 85, the Pedum Firm provides long lasting functional control. Heel raises and medial heel wedges can be used facilitating rapid clinical adjustment.



Due to it's unique material, we have designed adhesive backed heel raises (4mm & 6mm) and adhesive backed 4 degree heel wedges.

Code	Size
FIRM	From shoe size Jnr 12 - adult 11'

### Pedum Heel Raises

Each heel raise runs the full length of the calcaneum thus avoiding any 'rocking'



Size
4mm - Sizes 4.5 - 11
6mm - Size 4.5 - 11

### Pedum Medial Heel Wedges

4 degree medial rearfoot wedges. Designed to fit all Pedum insole patterns.



Size
Sizes 4.5 - 11

## MISCELLANEOUS INSTRUMENTS

## HOW TO TRADE INTERNATIONALLY WITH BAILEY INSTRUMENTS

### Heavy Duty Scissors



Code	Size
CH70	20cm

### Felt Scissors



Code	Size
CH129	20cm

### Tissue Forceps

(NIPPER ACTION)



Code	Size
CH04	10cm

### Tissue Forceps

(SCISSOR ACTION)



Code	Size
CH05	10cm

### Toenail Scissors



Code	Size
SD7-0271	11cm

### Reach your Feet

INNOVATIVE PRODUCT DESIGNED BY A PODIATRIST. THE LONG HANDLED 'GRIPPER' HOLDS THE TOE IN POSITION WHILST FILING WITH THE LONG-HANDLED FILE. 45CM LONG.



Code
CH146

### Contact Us

Export Sales Director - **Gemma Bailey**

T - **+44 7545 915 578**

E- [gemma@baileyinstruments.co.uk](mailto:gemma@baileyinstruments.co.uk)

### Customer Service

T - **+44 161 872 8707**

E- [sales@baileyinstruments.co.uk](mailto:sales@baileyinstruments.co.uk)


At Bailey Instruments, we operate our business with a heavy focus on Quality and Regulatory compliance which can be demonstrated through our product portfolio and customer service.

We are accredited to ISO13485:2016 with BSi Netherlands (CE 2797) and hold CE Certification on all our Medical Devices. We have successfully registered our products across the Globe, and work with many Importers and Distributors - Feel free to get in touch with our Export Director, Gemma Bailey, to learn more.

However, we don't believe compliance should just end there; therefore, we have strict procedures to ensure fair working conditions throughout our supply chain.

We strongly believe in improving lifestyle, whether this is our employees, suppliers & customers, health professionals, or patients!





Follow us on  

Unit 19 Guinness Rd Trading Estate,  
Guinness Rd, Trafford Park,  
Manchester, England, UK, M17 1SB

+44 (0)161 872 8707 | [sales@baileyinstruments.co.uk](mailto:sales@baileyinstruments.co.uk)  
 @BaileyInstLtd / @SUSOL  
[www.baileyinstruments.co.uk](http://www.baileyinstruments.co.uk) | [www.susol.co.uk](http://www.susol.co.uk)

**Bailey  
Instruments**  
investing in  
quality assurance



MD642415